



**MORE LIGHT**



## System requirements – JENOPTIK GRYPHAX<sup>®</sup> software & USB 3.0 microscope cameras

Beside of the latest software version on the memory card (included in delivery of the JENOPTIK GRYPHAX<sup>®</sup> cameras), we offer a free download after registration.

### JENOPTIK GRYPHAX<sup>®</sup> microscope camera software

**DOWNLOAD**

**Note 1:** For software and driver installation temporary administrator permission is necessary. Please contact your local IT department in advance!

**Note 2:** For software operation it is required permission of „full control“ for system folder & subfolders of OS:

Windows	"C:\ProgramDataJenoptik\..."
MacOS	"MacHD/Library/Application Support/Jenoptik\..."
Linux	"/var/lib/Jenoptik/..."

#### RECOMMENDED SYSTEM REQUIREMENTS (valid as of GRYPHAX version 2.3.0)

Personal Computer	Intel 11 <sup>th</sup> Gen i7 (min. 3.0 GHz) processor   16 GB RAM (dual-channel*) Apple Mac mini M4 2024   16 GB RAM
Operating System	Windows 10 version 22H2 or newer / Windows 11 22H2 – 64 bit MacOS Sequoia – 64 bit Linux Ubuntu 18.04 LTS – 64 bit
Data Interface	USB 3.0 integrated over PCI Express V2.0 (Renesas chip set)
Graphic Interface	Dedicated Graphic card equipped with on-board video memory
Monitor Resolutions	1920x1200 for AVIOR & SUBRA & KAPELLA & RIGEL & WEGA & POLARIS & PROKYON 3840x2160 for BETRIA & ALTAIR & ARKTUR & NAOS

\* dual-channel memory architecture mandatory to reach full camera performance

#### MINIMUM SYSTEM REQUIREMENTS

Personal Computer	Intel i5 (multi-core) processor or comparable   4 GB RAM (dual-channel*) Apple Mac mini M2 Pro 2023   8 GB RAM
Operating System	Windows 10 version 22H2 or newer / Windows 11 22H2 – 64 bit MacOS Sequoia – 64 bit Linux Ubuntu 18.04 LTS – 64 bit
Data Interface	USB 3.0 over PCI Express V1.1 (Renesas chip set)
Monitor Resolutions	1280x720 or higher

Release date: 2024-12-02

GRYPHAX<sup>®</sup> SERIES

JENOPTIK | Advanced Photonic Solutions  
JENOPTIK Optical Systems GmbH  
Goeschwitzer Strasse 25  
07745 Jena, Germany



**MORE LIGHT**

**NOT SUPPORTED:**

JENOPTIK GRYPHAX® cameras on USB2.0 interface! Connection issues & bad images can occur!  
JENOPTIK GRYPHAX® cameras under Windows 7 / 8 or older and 32 bit operating systems!

Learn more about JENOPTIK GRYPHAX® software and the easy workflows and tools.

Watch our video tutorials.

[VIDEOTUTORIALS](#)

HELPFUL? We appreciate your feedback.

[FEEDBACK](#)

We looking forward to hearing from you and remain sincerely,

Your JENOPTIK GRYPHAX® Team

Release date: 2024-12-02

GRYPHAX® SERIES

JENOPTIK | Advanced Photonic Solutions  
JENOPTIK Optical Systems GmbH  
Goeschwitzer Strasse 25  
07745 Jena, Germany