



Versatility –
giving you the freedom to work
with your favorite equipment.



JENOPTIK PROGRES GRYPHAX® & ProgRes® cameras

Hard- and software compatibility - overview.



JENOPTIK microscope cameras

hard- & software compatibility - overview

		PROGRES GRYPHAX® software			CapturePro® software	
		WIN	MAC	LINUX	WIN	MAC
USB 3.0	PROGRES GRYPHAX® SUBRA	X	X	X		
	PROGRES GRYPHAX® ARKTUR	X	X	X		
	PROGRES GRYPHAX® NAOS	X	X	X		
	PROGRES GRYPHAX® RIGEL*	X	X	X		
	PROGRES GRYPHAX® KAPELLA*	X	X	X		
	PROGRES GRYPHAX® PROKYON*	X	X	X		
USB 2.0	ProgRes® CT3 ^{USB}				X	
	ProgRes® CT5 ^{USB}				X	
	ProgRes® SpeedXT ^{core 3}				X	
	ProgRes® SpeedXT ^{core 5}	X			X	
	ProgRes® MF ^{USB}				X	
	ProgRes® CF ^{USB}				X	
FireWire 400	ProgRes® CT3				X	X
	ProgRes® C3				X	X
	ProgRes® C5				X	X
	ProgRes® C7				X	X
	ProgRes® MF				X	X
	ProgRes® CF				X	X
	ProgRes® CF ^{cool}	X	X	X	X	X
	ProgRes® MF ^{cool}	X	X	X	X	X
	ProgRes® CF ^{scan}				X	X
	ProgRes® MF ^{scan}				X	X
	ProgRes® C14 ^{Plus}	X	X	X	X	X

*camera model supported by GRYPHAX software version 1.1.4 or newer

If you need any further information please send an email to:

progres@jenoptik.com

Your PROGRES GRYPHAX® Team